




These products have been assessed by Arthritis Ireland's panel of patient experts as being Easy to Use for people with arthritis.

Please note that these products should only be used to treat the arthritis or condition as indicated on each.



KneeHab XP

Thigh muscle wastage can be a significant problem for arthritis sufferers as strengthening exercises are often restricted due to pain, swelling and lack of function. KneeHab XP works by strengthening and re-educating the thigh muscles using a unique, garment based, neuromuscular electrical stimulation device. Clinical research has proven that KneeHab increases muscle strength, decreases pain and aids in speeding your return to normal daily activities which has a positive impact on your quality of life. KneeHab's use pre Total Knee Replacement has shown that significant strength and functional gains can be achieved versus exercise alone and this could potentially delay surgery or if surgery is required

speed recovery post surgery.

Quattromed XP

Quattromed XP is a light-weight brace designed for arthritis patients suffering from medial knee pain (the medial side of the knee is the side closest to the other knee). It works by dissipating pressure above and below the knee to unload the knee joint thereby relieving pain, making it easier to walk and allowing you to enjoy a more active lifestyle.



Baxolve XP

Baxolve XP is a drug-free, non-invasive therapy which provides immediate and prolonged



back pain relief for arthritis sufferers. It combines two clinically proven therapies, TENS and Lumbar Support, in one integrated belt to support the spinal and abdominal muscles while at the same

time relieving pain and reducing the need for drug therapy. It is an easy to use home based therapy which can be worn during daily activities to help you regain mobility and quality of life.

Meditens XP

Meditens XP is a clinically proven, safe and effective means of providing a drug-free method of pain relief for arthritis sufferers. TENS fights pain by applying mild electrical pulses through pads placed on the skin blocking the pain signal to the brain and producing endorphins, the bodies own natural pain relievers. It is ideal for home use and has been proven to combat the pain experienced in arthritic joints allowing you to carry out your normal daily activities and improve your quality of life.





These products have been assessed by Arthritis Ireland's panel of patient experts as being Easy to Use for people with arthritis.

Please note that these products should only be used to treat the arthritis or condition as indicated on each.



Lumbofix

Lumbofix is a garment used in conjunction with Meditens XP or Medistim XP to isolate and treat lower back pain. It can be worn under clothing and ensures easy electrode placement at home to provide error free, repeatable treatment that results in better treatment outcomes for you.

Medistim XP

MediStim XP is an Electrical Muscle Stimulation device that re-educates and strengthens weakened muscles that may result from pain, swelling or reduced mobility. The device can be used to rehabilitate any muscle group thereby helping patients recover muscle strength, increase function and return to normal daily activities as soon as possible. It can also aid in the management of acute pain with its TENS program which is particularly effective in the treatment of arthritic pain. It is a home based therapy that offers patients comfort and control.



All products are available for sale from:

Freephone: 1800 511 511

www.neurotech.ie

info@neurotech.ie



Or log on to Arthritisireland.ie / online shop

